

HOW IT WORKS

HTS-flow is based on the IFOM/IEO/IIT Campus cluster, and therefore a Campus cluster account is needed to log in.

If you do not have a Campus cluster account yet, you have to ask on to the Service Desk system:

<http://servicedesk.ieo.it/helpdesk/WebObjects/Helpdesk.woa>

You can contact the service desk if you cannot remember your username or password.

Once you have obtained a cluster account, you need to contact the HTS-flow responsible person to be added in the system.

PRIMARY ANALYSIS

Submitting a job

From this page a computational biologist can run a primary analysis on a group of samples. Primary analyses consist of quality control of the raw reads, followed by filtering and alignment to reference genome.

The genomes currently available in HTS-flow are:

mm9 - Mus musculus

mm10 - Mus musculus

hg18 - Homo sapiens

hg19 - Homo sapiens

rn5 - Rattus norvegicus

dm6 - Drosophila melanogaster

Sequencing technologies currently implemented in HTS-flow are RNA-Seq, ChIP-Seq, DNaseI-Seq, BS-Seq.

SECONDARY ANALYSIS

Secondary Analyses are datatype-specific. HTS-flow supports:

Expression Quantification

Differential Genes Expression (DEG calling)

Peak Calling

Footprint Calling

This page is used to launch secondary analyses. In the first panel, the user can choose the appropriate type of secondary analysis, while in the second panel the available aligned samples can be chosen by clicking on them and then pressing the “SELECT” button. A final panel, specific for the chosen secondary analysis, will open automatically.

When a secondary analysis is complete, it will be shown in the COMPLETED ANALYSIS page, under the SECONDARY JOBS table. The output of each analysis is stored in the specific folder within HTS-flow. Secondary analysis results can then be accessed in two ways:

(1) With a web browser, at the URL:

http://www.bioinfo.ieu.eu/BAgroup/HTS-flow/DB/secondary/SECONDARY_ID/

where SECONDARY_ID corresponds to the Secondary ID assigned by HTS-flow.

(2) On the campus cluster (grid.ieu.eu) at the address:

/data/BA/public_html/HTS-flow/DB/secondary/SECONDARY_ID

where SECONDARY_ID corresponds to the Secondary ID assigned by HTS-flow.

The output of a secondary analysis is saved in RDS format, which can be loaded within R with the function readRDS:

```
results=readRDS('path_to_your_result_file')
```

Refer to <https://stat.ethz.ch/R-manual/R-devel/library/base/html/readRDS.html>

for more information about this function.

Expression Quantification

Submitting jobs

After choosing the samples, the user must fill a table with two fields:

SAMPLE: the PRIMARY ID of the sample;

MIX: the ERCC spike-in Mix (either 1 or 2) added for normalization. If no mix was added, leave empty. This feature has not yet implemented.

The ADD and REMOVE buttons can be used to create/delete more lines in the table and quantify the expression of multiple samples within the same secondary analysis.

Finally, the job can be submitted by clicking the SUBMIT button (Figure 1). HTS-flow will assign a SECONDARY ID to the job and automatically redirect the web browser to the RUNNING ANALYSES page, which contains the list of running jobs. When the analysis is complete, it will be moved to the COMPLETED ANALYSES webpage, in the SECONDARY JOBS panel.

EXPRESSION QUANTIFICATION

SELECTED SAMPLES

PRIMARY ID	SAMPLE ID	sample_name	READS NUM	OPTIONS	ref_genome_atn	METHOD	SOURCE	USER
2120	81	Sample_S_31BD_2h_EtOH_S7704	32659645	view/hide	mm9	RNA-Seq	EXTERNAL	vbianchi
2119	82	Sample_S_31BD_2h_100nM_OHT_S7705	30594802	view/hide	mm9	RNA-Seq	EXTERNAL	vbianchi
2118	83	Sample_S_45355_2h_EtOH_S7706	26959743	view/hide	mm9	RNA-Seq	EXTERNAL	vbianchi
2117	84	Sample_S_45355_2h_100nM_OHT_S7707	23825746	view/hide	mm9	RNA-Seq	EXTERNAL	vbianchi

For Expression Quantification you have to fill the following form with IDs in order to have:

- in SAMPLE box the ID of the sample;
- in MIX box the type of ERCC Mix used if available;

SAMPLE	MIX	ADD	REMOVE
2120			
2119			
2118			
2117			

If you feel confident you can submit this job!

SUBMIT

Figure 1: HTS-flow Expression Quantification (Secondary Analysis) web interface. In this example, 4 samples have been selected for expression quantification analysis. The primary IDs were used for filling the table, while no spike-in Mix was selected.

Output

The Expression Quantification analysis will create two different output files: RPKMS.rds and eRPKMS.rds.

(1) RPKMS.rds: absolute quantification of gene expression in terms of Reads per Kilobase per Million of mapped reads (RPKM).

(2) eRPKMS.rds: absolute quantification of gene expression in terms of Reads per Kilobase per Million of mapped exonic reads (eRPKM).

Both files are R data frames where each row is a gene (row names are Gene Symbols) and each column is a sample, and values correspond to gene expressions. In Figure 2, the first 6 rows of the RPKMS.rds file from a two-sample analysis are shown; the column names are the sample names retrieved by the LIMS or defined by users when submitting external data to HTS-flow.

	p53K0_1_4h_7Gy_S_B220plus	p53K0_2_4h_7Gy_S_B220plus
0610007P14Rik	10.7384306	7.59893422
0610009B22Rik	2.3839230	2.33540480
0610009D07Rik	23.2899388	14.05794407
0610009020Rik	3.5231083	2.81864757
0610010B08Rik	0.2428211	0.09795024
0610010F05Rik	1.4856537	1.16142305

Figure 2: HTS-flow Expression Quantification output. The first 6 rows of the result table are shown. For each gene, the RPKM value is reported in each sample.

Differential Gene Expression (DEG calling)

Submitting jobs

Calling Differentially Expressed Genes (DEGs) is performed with DESeq2. The experimental design requires two conditions, typically a treated set of samples and a control set of samples. The presence of replicates in at least one condition is essential.

After choosing the samples, the user must fill a table with three fields:

SAMPLE: the PRIMARY ID of the sample.

CONDITION: one of the two classes of samples (e.g.: treated and control) in the experimental design.

MIX: the ERCC spike-in Mix (either 1 or 2) added for normalization. If no mix was added, leave empty. This feature has not yet implemented.

EXP NAME: a name for the analysis, to be used for the output file. Avoid spaces, use instead underscores.

The ADD and REMOVE buttons can be used to create/delete more lines in the table and add replicates to the conditions.

Finally, the job can be submitted by clicking the SUBMIT button (Figure 3). HTS-flow will assign a SECONDARY ID to the job and automatically redirect the web browser to the RUNNING ANALYSES webpage, showing the list of running jobs. When the analysis is complete, the analysis will be visible on the COMPLETED ANALYSES webpage, under the SECONDARY JOBS panel.

For DEG analysis you have to fill the following form with IDs in order to have:

- in SAMPLE box the ID of the sample;
 - in CONDITION box provide the name for the condition (usually **treat** and **control**)
- Please bear in mind that the names are used by DESeq2 in alphabetical order. So if you label two conditions 'a' and 'b' the analysis will be performed 'b vs a';
- in MIX box the type of ERCC Mix used if available (is yet in testing);
 - do not use space when you name analysis or labels, use instead the underscore (_).
 - do not use numbers as starting character for naming the analysis or labels (ex. 0h_Myc, use instead Myc_0h).

SAMPLE	CONDITION	MIX	ADD	REMOVE	EXP NAME
2202	ctrl				mdr2_ko_TvsC
2203	ctrl				
2204	ctrl				
2205	ctrl				
2206	ctrl				
2212	tumor				
2213	tumor				
2214	tumor				
2215	tumor				
2216	tumor				

If you feel confident you can submit this job!

SUBMIT

Figure 3: HTS-flow DEG web interface. The experiment has two conditions: ctrl and tumor, which were associated to the primary IDs of the samples. An experiment name was provided (mdr2_ko_TvsC).

Output

Differential Gene Expression analysis will create a single output file, corresponding to the EXP NAME field, in the form 'EXP NAME'.rds.

This file contains an R data frame where each row is a gene (Gene Symbol) and columns list the DESeq2 default outputs (see Figure 4):

baseMean: the mean gene expression over all samples in the two conditions

log2FoldChange: log2 Fold Change, treated vs untreated

lfcSE: standard error, treated vs untreated

stat: Wald test statistic

pvalue: Wald test p-value

padj: Benjamini-Hochberg adjusted p-values (False Discovery Rate)

	baseMean	log2FoldChange	lfcSE	stat	pvalue
0610007P14Rik	332.1377	-0.05977906	0.13356388	-0.4475691	0.65446424
0610009B22Rik	197.6424	-0.01904953	0.16860606	-0.1129825	0.91004445
0610009020Rik	790.7445	-0.12799817	0.08295034	-1.5430700	0.12281379
0610010B08Rik	278.4833	0.38895479	0.15269362	2.5472890	0.01085635
0610010F05Rik	853.9069	0.21065296	0.09214272	2.2861595	0.02224493
0610010K14Rik	205.6984	-0.21230918	0.10988940	-1.9320260	0.05335630
0610007P14Rik	0.8806745				
0610009B22Rik	0.9715694				
0610009020Rik	0.4764572				
0610010B08Rik	0.1659328				
0610010F05Rik	0.2296517				
0610010K14Rik	0.3410278				

Figure 4: DEG calling output in HTS-flow (first 6 rows). Rows correspond to genes, while the different DESeq2 outputs are on the columns.

Peak Calling

Submitting jobs

Peak Calling is performed with MACS2.

After choosing the samples, the user must fill a table with three fields:

SAMPLE1: the PRIMARY ID of the sample used as input in the peak call.

SAMPLE2: the PRIMARY ID of the sample where peaks should be called.

LABEL: label associated to SAMPLE2, which will be used as reference name for SAMPLE2 in this secondary analysis. Avoid spaces, use instead underscores.

The ADD and REMOVE buttons can be used to create/delete more lines in the table and performing more peak calls within the same secondary job.

A name for the analysis must to be assigned through the **EXP NAME** input form. Avoid spaces in the name, use instead underscores.

The last part of the form allows to select the parameters for the peak call performed by MACS2.

Peak shapes must chosen accordingly to the specific ChIP-seq experiment: for example, NARROW peaks for transcription factors, BROAD peaks for histone marks.

With the NARROW/BROAD option both calls will be performed and the union of peaks from the NARROW and BROAD analysis will be output. This option will affect the annotation analysis performed by HTS-flow on the identified peaks as explained below.

Finally, the job can be submitted by clicking the SUBMIT button (Figure 5). HTS-flow will assign a SECONDARY ID to the job and automatically redirect the web browser to the RUNNING ANALYSES webpage, showing the list of running jobs. When the analysis is complete, the analysis will be visible on the COMPLETED ANALYSES webpage, under the SECONDARY JOBS panel.

For peak calling analysis you have to fill the following form with IDs in order to have:

- in SAMPLE 1 box the ID of the input;
- in SAMPLE 2 box the ID of the ChIP;
- if you do not have an input for the ChIP fill both SAMPLE 1 and SAMPLE 2 with the same ID.
- if you have more than one ChIP click the ADD button to insert another ChIP.
- do not use space when you name analysis or labels, use instead the underscore (_).
- do not use numbers as starting character for naming the analysis or labels (ex. 0h_Myc, use instead Myc_0h).

All fields are mandatory.

SAMPLE 1	SAMPLE 2	LABEL	ADD	REMOVE
2122	2224	Anti_dLsd_1		
2123	2225	Anti_dLsd_2		

EXP NAME

YOU WANT TO FIND

P-value	<input type="text" value="0.00001"/>	e.g.: 0.00001 is 10e-5
Options	<input type="text" value="-mfold=7,30"/>	

If you feel confident you can submit this job!

Figure 5: HTS-flow Peak Calling web interface.

Output

The output of the peak calling analysis is distributed in two folders: the NARROW/ or BROAD/ folders (depending on the type of call) contain the MACS2 output, i.e. a bed file containing the genomic locations for each peak identified for each sample and a saturation table file for each sample. Saturation reports at different fold-enrichments, the proportion of peaks that could still be detected when using 80% to 20% of the sequence reads. Saturation file name is in the form 'LABEL'_saturation.txt.

Besides, in the annotation/ folder the Peak Calling analysis will create an output file per sample submitted, whose name will be in the form 'LABEL.rds'. Each file contains a GRanges object where each element is a genomic interval (a peak), complemented by information obtained with the GRannotate and GRenrichment functions from compEpiTools

A SECONDARY JOBS

454	peak_calling	2015-05-12 12:01:04	2015-05-14 14:09:16	50:08:12	/data/BA/public_html/HTS-flow/DB/secondary/454/	Link	view/hide	EuMyC DNase I high-depth DHSs against merged input	pbora
453	peak_calling	2015-05-12 11:09:30	2015-05-12 16:45:12	06:55:36	/data/BA/public_html/HTS-flow	Link	view/hide		

B Index of /BAgroup/HTS-flow/DB/secondary/454/NARROW

Name	Last modified	Size	Description
Parent Directory	-	-	-
EuMyC_C_DHS_20_model.r	13-May-2015 13:36	56K	
EuMyC_C_DHS_20_peaks.bb	13-May-2015 15:18	699K	
EuMyC_C_DHS_20_peaks.bed	13-May-2015 15:18	1.7M	
EuMyC_C_DHS_20_peaks.encodePeak	13-May-2015 15:18	2.5M	
EuMyC_C_DHS_20_peaks.xls	13-May-2015 15:18	2.3M	
EuMyC_C_DHS_20_pg_table.txt	13-May-2015 14:58	582K	

C Index of /BAgroup/HTS-flow/DB/secondary/454/annotation

Name	Last modified	Size	Description
Parent Directory	-	-	-
EuMyC_C_DHS.rds	14-May-2015 11:11	1.7M	
EuMyC_P_DHS.rds	14-May-2015 12:24	1.7M	
EuMyC_T_DHS.rds	14-May-2015 12:40	1.9M	

Figure 6: Peak calling output in HTS-flows. (A) Output for secondary ID 454 in the completed analyses page. Clicking on the "Link" button opens the folders with the results: MACS2 output (B) and annotated rds files (C).

For each peak, the following fields are available:

- enrichment:** $\log_2(\text{ChIP}/\text{N1-input}/\text{N2})$, where ChIP is the number of reads falling in the interval in the sample, N1 is the library size of the sample, input is the number of reads falling in the interval in the input and N2 is the library size of the input. Computed with the GRenrichment function from compEpiTools.
- summit:** position of maximum coverage of the peak. Computed with the GRcoverageSummit function from compEpiTools.
- midpoint:** the midpoint of the peak. Computed with the GRmidpoint function from compEpiTools.

Annotation is computed in two distinct ways, depending on the type of peak calling requested:

- NARROW calls computes annotation from the summit of the peak
- BROAD calls computes annotation from the midpoint of the peak.
- NARROW/BROAD calls computes annotation from the midpoint of the peak.

Output

Wellington outputs a bed file containing the genomic locations for each footprint for each sample. These file are located in the output folders footprints/. The bed file contains the the genomic locations for each genomic footprint, followed by a score associated by the footprint caller assessing its statistical significance.

Besides reporting the original wellington output (a txt file), HTS-flow converts the bed files in R GRanges objects in the form 'EXP_NAME'.rds in the output folder footprints/.

A SECONDARY JOBS

Foot										Foot					
SECONDARY ID	METHOD	START	END	TIME (hh:mm:ss)	FOLDER OUTPUT	LINK OUTPUT	INFORMATIONS				DESCRIPTION	USER			
460	footprint	2015-05-15 11:27:38	2015-05-15 15:52:50	04:25:12	/data/BAGpublic_html/HTS-flow/DB/secondary/460/	Link	<div>view/hide</div>						mycER 3T9 footprints merged input Edit	pbora	
							FOOTPRINT ID	SECONDARY ID	EXP NAME	REF SEC ID	CALLER	PVALUE			DESCRIPTION
							19	460	mycER 3T9 footprints	450	wellington	0.0000000001			mycER 3T9 footprints merged input
							20	460	mycER 3T9 footprints	450	wellington	0.0000000001			mycER 3T9 footprints merged input
451	footprint	2015-05-12 11:16:23	2015-05-13 13:54:33	26:38:10	/data/BAGpublic_html/HTS-flow/DB/secondary/451/	Link	<div>view/hide</div>						MycER 3T9 DNase I footprints 0h and 4h Edit	pbora	
449	footprint	2015-05-11 11:24:04	2015-05-11 11:26:02	00:01:58	/data/BAGpublic_html/HTS-flow/DB/secondary/449/	Link	<div>view/hide</div>						EuMyc C.P footprints Edit	pbora	

B Index of /BAGroup/HTS-flow/DB/secondary/460/footprints

Name	Last modified	Size	Description
Parent Directory			
DNaseI_0h_footprints.rds	14-May-2015 22:57	4.2M	
DNaseI_0h_footprints.txt	14-May-2015 22:57	19M	
DNaseI_4h_footprints.rds	15-May-2015 15:27	4.0M	
DNaseI_4h_footprints.txt	15-May-2015 15:26	19M	

Apache/2.2.22 (Debian) Server at www.bioinfo.leo.eu Port 80

C

chr1	3585829	3585855	Unnamed128627	-42.24478149656517	+
chr1	3585100	3585116	Unnamed128628	-13.108855884707419	+
chr1	4768322	4768344	Unnamed128629	-19.10218744293243	+
chr1	4775485	4775527	Unnamed128630	-27.28617814152878	+
chr1	4775567	4775579	Unnamed128631	-13.789649756271892	+
chr1	4775627	4775664	Unnamed128632	-26.91556636820937	+
chr1	4775668	4775694	Unnamed128633	-18.994383457289984	+
chr1	4775728	4775769	Unnamed128634	-53.569175879680846	+
chr1	4775801	4775875	Unnamed128635	-148.36808164623905	+
chr1	4775912	4776088	Unnamed128636	-25.26126814706762	+
chr1	4776041	4776057	Unnamed128637	-10.637658811816426	+
chr1	4776081	4776093	Unnamed128638	-18.55395114315376	+
chr1	4797736	4797752	Unnamed128639	-18.78633186648126	+
chr1	4797814	4797834	Unnamed128640	-19.426191375537957	+
chr1	4797876	4797916	Unnamed128641	-29.21439248924102	+
chr1	4797963	4797977	Unnamed128642	-13.036669414998667	+
chr1	4798015	4798065	Unnamed128643	-20.29775592638477	+

Figure 8: Footprint results on HTS-flow. (A) Output for footprint analysis with secondary ID 460 in the completed analyses page. Clicking on the "Link" button opens the folders with the results: bed file with footprints identified from the analysis is available both in rds and txt format. In (C) is shown the header of the bed file containing footprints genomic locations.

MERGING ANALYSIS

In this page a user can select a group of aligned samples and pool their reads to obtain a merged alignment file. To be merged, samples need to be aligned to the same reference genome.

Querying completed jobs

The top table (MERGED SAMPLES) shows the completed merging jobs. Each line of this table consists in a sample used in a merging analysis. To identify all the samples used in a merging job, gather together all the rows with the same 'MERGE ID' (for example, by filling the desired 'MERGE ID' on the top of the table).

Submitting a job

The bottom panel (MERGING) can be used to select samples to be merged. Clicking on the SELECT button opens a new panel with two fields:

SAMPLE NAME: the name for the merged sample.

REMOVAL OF DUPLICATES: default TRUE: removal of duplicates in the merged file (suspected to be PCR duplicates). For DNaseI-Seq at high depth (footprint calling), change it to FALSE.

Finally, the job can be submitted by clicking the SUBMIT button. HTS-flow will assign a MERGE ID to the job and automatically redirect the web browser to the RUNNING ANALYSES webpage, showing the list of running jobs. When the analysis is complete, the analysis will be visible on the COMPLETED ANALYSES webpage, under the SECONDARY JOBS panel.